Louisiana Office of Public Health Laboratories	
Test Name	E. coli Enumeration Test of Source Water Samples
PHL Location	Amite Laboratory, Central Laboratory
CPT Code	N/A
Synonyms	LT2
Brief Description of Test	The enumeration of <i>E. coli</i> in source water using SM9223 (Colilert Quanti-tray 2000).
Possible Results	E. coli Most Probable Number (MPN)/100 ml
Reference Range	A mean of 50 <i>E. coli</i> for flowing stream and a mean of 10 <i>E. coli</i> for reservoir/lakes.
Specimen Type	Source Water
Specimen Container(s):	Sterile, disposable 7oz. Nasco Whirl-Pak bag. Only containers QCed and supplied by the laboratory are suitable for use. DO NOT purchase privately.
Minimum volume accepted:	120 ml
Collection Instructions	 Collect representative sample according to established protocols. Use a Whirl-Pak bag and fill leaving at least 1 inch of air space at top of container. Whirl bag to close. Do not twist tabs together.
Storage and Transport Instructions	Store & transport sample(s) in clean iced cooler(s). Samples must be at <10°C when received. Maximum holding time is 30 hours.
Causes for Rejection	 Sample >30 hours old. Frozen water sample. No sample ID. Unsatisfactory collection container. Insufficient amount of sample. Sample container leaking. Sample information form not complete. Sample temperature <10°C upon receipt.

Limitations of the Procedure	Not suitable for marine water or any sample altered by any pre- enrichment or concentration.
Interfering Substances	Samples containing humic material with innate color.
References	 Idexx Laboratories Colilert Colilert 18 package inserts LA DHH-OPH Laboratory Services SOP DOLSOPDWM-01 Chromogenic Fluorogenic Substrate Method.
Additional Information	Complete a Lab 8 Laboratory Request and Report Form. The following is required Lab 8 information: Sample sequence #s Project Code Name of Supply Collected by (name) PWS ID Date/Time Collected Type of source water Turbidity
Release Date	11/2013

Warning: If you have printed a copy of this information please be advised that the Louisiana Office of Public Health Laboratories website and methods are updated on a regular basis. Please check the on-line version of this document to ensure you are relying on the most recent release.

LO.FM.GEN.043 V2 4 2013